

CAVS® Cart description

Mono-block chassis manufactured in HDPE. High rigid structure resistance, designed to lodge the maximum equipment in a minimum space (570 x 430 x 930 mm). The own structure has different features: ergonomic handles, upper tray for equipment, lodgings for dripping support, space specifically designed for documents or a folding tray, lodging for an O₂ cylinder of 5 litres (Ø140), standard supports for containers, 4 non-skid wheels of 100mm (1 of them antistatic and 2 others coming with brakes).

Modular chest with removable drawers, mobile dividers and adhesive labels, protected by a Cardiac Massage Board in transparent methacrylate, including a plastic breakaway seal. It serves as protection against uncontrolled openings of drawers. The content of the drawers (medicines and components) is not provided with the CAVS (see the Listing of Recommended Equipment in the CAVS instructions manual).

Other optional equipments	
Electrical Suction Units	Ventilators
<div></div> <div>5320236 V7mx (AC+DC+battery)</div> <div>5320261 V7plusb (DC+battery)</div>	<div></div> <div>5200800 Vitae 40</div>
<div></div> <div>4320033 AC/DC charger (for V7 plusb)</div>	

Cart for Advanced
Vital Support

CAVS®

The most compact,
arranged and functional system
in the market



The CAVS concept ...making emergency carts since 1985




Basic Unit

5230500

Mono-block chassis with 4 wheels and 7 drawers: 3 of 350 x 150 x 100 mm and 4 of 350 x 150 x 75 mm	
Resuscitator container:	4230502
• Adult resuscitator Revivator Plus (silicone + PSU) with face mask no. 5 and O ₂ reservoir 2,5 L	5200601
• Child resuscitator Revivator Plus (silicone + PSU) with face mask no. 3	5200612
Collection jar for suction (PC 1 L)	4383001
EN support for waste recipient	4383015
Waste recipient (yellow 1 L)	3230515
Cardiac massage board (Works also as a protection against uncontrolled openings of drawers)	3230505
Dripping support	4230506
File - Folder for documentation	3230507
Instructions manual	9230500
Adhesive label kit for drawers	3230508
Plastic breakaway seal (x5)	3230510

O₂ + Suction (optional)



Kits (optional)



CAVS cart as shown in the cover is composed by the Basic unit and O₂ + Suction 5230500 + 5302011 + 5505011 + 0510003