

Laringoscopios

Laryngoscopes



Palas con Fibra Óptica Integrada

Blades with Fiber Optic Internal light conductor

Lames avec Fibre Optique intégrée



Palas de acero inoxidable, AISI-304, no magnéticas, con fibra óptica integrada.

Stainless steel blades, AISI-304, non-magnetic, with fiber optic internal light conductor.

Lames en acier inoxydable, AISI-304, non magnétiques, avec fibre optique intégrée.

Megaray

		a)	b)
0130120	Macintosh # 0	81 mm	55 mm
0130121	Macintosh # 1	94 mm	70 mm
0130122	Macintosh # 2	16 mm	90 mm
0130123	Macintosh # 3	135 mm	115 mm
0130124	Macintosh # 4	155 mm	132 mm

0130125 Estuche rígido con: mango Standard, palas # 1-2-3-4 y bombilla de repuesto.

Rigid case with: Standard handle, blades # 1-2-3-4 and spare bulb.

Mallette rigide avec: manche Standard, lames # 1-2-3-4 et ampoule de rechange.

Flexiray



Mayor flexibilidad y control para INTUBACIÓN DIFÍCIL.


Greater flexibility and improved control for DIFFICULT CASES

Plus grande flexibilité et contrôle pour INTUBATION DIFFICILE

		a)	b)
0130129	Macintosh # 3	135 mm	115 mm
0130130	Macintosh # 4	155 mm	132 mm



0130128 Estuche rígido (vacío) para 1 mango y 2 palas.
Rigid case (empty) for 1 handle and 2 blades.
Mallette rigide (vide) pour 1 manche et 2 lames.





Palas de acero inoxidable, AISI-304, no magnéticas, con fibra óptica extraíble.

Stainless steel blades, AISI-304, non-magnetic, with removable fiber optic light.

Lames en acier inoxydable, AISI-304, non magnétiques, avec fibre optique détachable.

Palas con Fibra Óptica Extraíble

Blades with Removal Fiber Optic light

Lames avec Fibre Optique détachable

		a)	b)
0130110	Macintosh # 0	80 mm	55 mm
0130111	Macintosh # 1	92 mm	68 mm
0130112	Macintosh # 2	98 mm	75 mm
0130113	Macintosh # 3	30 mm	115 mm
0130114	Macintosh # 4	155 mm	130 mm
0130109	Miller # 00	68 mm	45 mm
0130105	Miller # 0	78 mm	55 mm
0130101	Miller # 1	102 mm	78 mm
0130102	Miller # 2	155 mm	130 mm
0130103	Miller # 3	195 mm	172 mm
0130104	Miller # 4	205 mm	182 mm

- 0130115** Estuche rígido con: mango Standard, palas Macintosh # 1-2-3-4 y bombilla de repuesto.
Rigid case with: Standard handle, Macintosh blades # 1-2-3-4 and spare bulb.
Mallette rigide avec: manche Standard, lames Macintosh # 1-2-3-4 et ampoule de rechange.

Mangos y accesorios F.O.

Handles and accessories F.O.

Manches et accessoires F.O.

Mangos Hook-on ISO 7376/3 "Green" System German, con bombilla luz de Xenon (latón cromado)

Hook-on ISO 7376/3 "Green" System German handles, with xenon light bulb (chromium plated brass)

Manches Hook-on ISO 7376/3 "Green" System German, avec ampoule lumière de xenon (laiton chromé)



0130117 Standard
Ø 28 mm x 160 mm



0130126 Penlight
Ø 18 mm x 160 mm

(para uso pediátrico)
(for pediatric use)
(pour usage pédiatrique)



0130127 Stubby
Ø 34 mm x 125 mm

(para obeso mórbido)
(for morbid obese)
(pour obèse morbide)

Bombilla Xenon 2,5 V • Xenon bulb 2,5 V • Ampoule Xenon 2,5 V

Estuche rígido • Rigid case • Mallette rigide

Estuche con cremallera • Pouch with zip • Trousse à fermeture éclair



0130116

0130108

0130010



(para 1 mango y 4 palas)
(for 1 handle and 4 blades)
(pour 1 manche et 4 lames)

(para 1 mango y 3 palas)
(for 1 handle and 3 blades)
(pour 1 manche et 3 lames)

Palas Convencional
Conventional Blades
Lames Conventiennel



Palas de acero inoxidable, AISI-304, no magnéticas, con bombilla de luz convencional.

Stainless steel blades, AISI-304, non-magnetic, with conventional light bulb.

Lames en acier inoxydable, AISI-304, non magnétiques, avec ampoule de lumière conventionnelle.

			a)	b)
0130003	Macintosh # 0	82 mm	60 mm	
0130004	Macintosh # 1	97 mm	75 mm	
0130005	Macintosh # 2	117 mm	95 mm	
0130006	Macintosh # 3	137 mm	115 mm	
0130007	Macintosh # 4	157 mm	135 mm	
0130008	Macintosh # 5	177 mm	155 mm	
0130025	Miller # 0	77 mm	55 mm	
0130026	Miller # 1	102 mm	80 mm	

0130001 Mango Standard • Standard handle • Manche Standard Ø 28 mm

0130002 Mango Penlight • Standard Penlight • Manche Penlight Ø 19 mm

0130012 Bombilla 2,5 V • Light bulb 2,5 V • Ampoule 2,5 V

0130010 Estuche cremallera • Pouch with zip • Trousse fermeture éclair

0130000 Estuche cremallera con mango Standard y palas Macintosh # 2-3-4
 Pouch with zip with Standard handle and Macintosh blades # 2-3-4
 Trousse avec manche Standard et lames Macintosh # 2-3-4



0130106 Mango Standard de plástico Hook-on
 ISO 7376/3 "Green" System German Ø 28 mm
 Standard plastic handle Hook-on
 ISO 7376/3 "Green" System German Ø 28 mm
 Manche plastique Standard Hook-on
 ISO 7376/3 "Green" System German Ø 28 mm

Palas metálicas DESECHABLES, de Fibra Óptica

Fiber Optic DISPOSABLE metallic Blades

Lames métalliques USAGE UNIQUE, Fibre Optique

			a)	b)
0130141	Macintosh # 0	88 mm	60 mm	
0130142	Macintosh # 1	92 mm	64 mm	
0130143	Macintosh # 2	98 mm	72 mm	
0130144	Macintosh # 3	130 mm	105 mm	
0130145	Macintosh # 4	155 mm	32 mm	



HERSILL, S.L.

Equipos Médicos - Medical Devices

Puerto de Navacerrada, 3 • Polígono Industrial Las Nieves
 28935 Móstoles, (Madrid), SPAIN
 Tel.: +34 91 616 41 11 • Fax: +34 91 616 48 92
 E-mail: info@hersill.com
 export@hersill.com

www.hersill.com



EN ISO 13485:2003



EN ISO 9001:2000



Nº de folleto
 Leaflet No.
 Brochure n°
 0130/08-09