

R-8000

Basic cart for Cardiopulmonary
Resuscitation



R-8002



HERSILL
Medical Devices

The **R-8000** offers you what you need: a cart that is big enough but not excessive, a central tray for monitors, very resistant, a basket under the tray and excellent castors for a noiseless, speedy and very

stable motion. It can be equipped with everything essential for RCP. However, its design allows for great versatility and, therefore, can be useful for many other functions, as an advanced life support trolley.

Basic composition

- Cart made in steel, painted with epoxy resin and equipped with 4 swivel castors of 125 mm, of which 1 is antistatic and 2 come with brakes.
- A module of 310 x 400 x 470 mm with two drawers (R-8000) and four drawers (R-8002), with dividers and key.
- Metallic basket over the drawers (only R-8000).
- Tray 450 x 400 mm, suitable to support up to 30 kg.
- Two side rails.
- Oxygen and suction equipment composed by: oxygen regulator, flowmeter 0-15 L/min, humidifier, venturi VacuJet suction device (50 cmHg – 650 mbar) and collection jar 1 litre.
- Support for two cylinders with fasteners.
- Manual resuscitator Revivator Plus for adults with two facemasks (adult and child).
- One IV Pole which can be placed at both sides of the cart.
- Dimensions: 1400 x 470 x 680 mm.
- Weight: 41 kg (R-8000) – 56 kg (R-8002).

R-8000

5230000

R-8002

5230100

Optional accessories

0510005 Oxygen cylinder 1400 litres.

0200020 Cardiac massage board.

5425900 Stainless steel probe support (Ø10 x 20 cm).

5200612 Manual resuscitator Revivator Plus child with face mask size 3.



5220100 R-2100 intubation set composed by:

- Nylon transport bag
- Laryngoscope: handle and 2 Macintosh blades size 2 and 4
- 6 Endotracheal tubes with balloon ID (mm) 5,6,7,8,9,10
- 2 Stylets for endotracheal tubes
- Serynge 100 ml for inflating tubes
- 3 Guedell tubes of different sizes
- Exploration lantern
- Cloth sticking plaster roll
- Magill forceps
- Pean forceps
- Bandage scissors



0300015 Manual suction pump.



5300200 VacUp suction foot pedal.



5320236 V7mx portable suction unit, 100/240 V_{ac} – 12 V_{dc} - internal battery.



5200800 Vitae 40 transport ventilator.



HERSILL

Medical Devices

HERSILL, S.L.

Puerto de Navacerrada, 3 • Polígono Industrial Las Nieves
28935 Móstoles (Madrid) SPAIN
Tel.: +34 91 616 60 00 • Fax: +34 91 616 48 92
info@hersill.com • export@hersill.com
www.hersill.com



ISO 13485



ISO 9001:2008



Leaflet No.
IN-523/14/07